

# DG 30



Research Microscopy for Professionals



LABO®

Worlds Within™

## SPECIFICATIONS



OPTICAL SYSTEM	Infinity Optical System (IOS)
STAND	C shaped, corrosion resistant, Anti mould treated, single mold modular stand
NOSEPIECE	Quintuple, Reversed (Inward Tilt) nosepiece, with rubber grip
OBSERVATION TUBE	30° Inclined Siedentopf Binocular tube, 360° rotatable with IPD 55-75 mm and dioptic adjustment
EYEPIECES	WF 10X/22mm with Eye-guards. Anti Mould Treated
OBJECTIVES	IP series Infinity corrected Plan 4x, 10x, 40x (SL), 100x (SL, Oil) , Anti Mould Treated
ILLUMINATION	Koehler Illumination, 6V 20W Halogen Lamp, Universal Power Supply 100 V to 240 V, AC 50-60 Hertz, CE Compliant
IMAGE DEVICE	3.0 Mega pixel Color CMOS sensor
RESOLUTION	2048(H) x 1536 (V) pixels, 20 fps at 1024x768 pixels
LIGHT SENSITIVITY	1 Lux
INTERFACE	USB ver 2.0
SOFTWARE	dgsoft ProMED ver 5.06 Image Analysis Software
SYSTEM REQUIREMENTS	CPU : Pentium IV or Higher, RAM : 512 MB or more, Hard Disk Space : 100 MB or more, Windows 98SE/ME/2000/XP, CD-Rom drive and USB 2.0 port



[www.labomicrosystems.com](http://www.labomicrosystems.com)

USA Facility:

### LABO Microsystems

6960-A Camino Maquiladora  
San Diego, CA 92154, USA.

Tel : +1 604 944 0622, Fax: +1 604 608 3225

Tel Toll Free USA/Canada 1 877 650 0117

E-mail: [sales@labomicrosystems.us](mailto:sales@labomicrosystems.us)

German Facility:

### LABO Microsystems GmbH

Bäckerstrasse 5  
21244 Buchholz, Germany

Tel: + 49 4181 9906360

Fax: + 49 4181 9906369

E-mail: [hq@labomicrosystems.com](mailto:hq@labomicrosystems.com)

Global Sales:

### LABO Microsystems

P.O. Box 213671, Jumeirah, Dubai.

United Arab Emirates.

Tel: + 97 150 242 0048

Tel: + 91 941 666 6666

E-mail: [labo@labomicrosystems.com](mailto:labo@labomicrosystems.com)



JAS-ANZ



ISO 9001:2000  
File No. 15292.04

