

DG Lab



Research Microscopy for Professionals



LABO®

Worlds Within™

SPECIFICATIONS



OBSERVATION TUBE	Trinocular 30° inclined, 360° rotatable, IPD 55-75 mm, with built-in camera and digitizer
EYEPIECES	WF 10x/18mm with Eye-guards, Anti Mould Treated. Optional WF 15x/16mm, WF 20x/12mm
OBJECTIVES	XP series Semi Plan Achromatic, Anti Mould Treated, 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65 SL, 100x/NA 1.25 SL, Oil. Optional 20x/ NA 0.40 SL, 60x/NA 0.85 SL
ILLUMINATION	Koehler Illumination, 6V 20W Halogen Lamp, with Universal Power Supply 100 V to 240 V, AC 50-60 Hertz, CE Compliant. Optional 6V 30W Halogen Lamp
IMAGE DEVICE	CCD Camera
RESOLUTION	Live/Still image resolution 720 x 576 pixels (PAL), 720 x 480 (NTSC)
LIGHT SENSITIVITY	1 Lux
INTERFACE	USB ver 2.0, S-Video & RCA
SOFTWARE	dgsoft ProMED ver 5.06 Image Analysis Software
SYSTEM REQUIREMENTS	CPU : Pentium IV or Higher, RAM : 512 MB or more, Hard Disk Space : 100 MB or more, Windows 98SE/ME/2000/XP, CD-Rom drive and USB 2.0 port



www.labomicrosystems.com

USA Facility:
LABO Microsystems
6960-A Camino Maquiladora
San Diego, CA 92154, USA.
Tel : +1 604 944 0622, Fax: +1 604 608 3225
Tel Toll Free USA/Canada 1 877 650 0117
E-mail: sales@labomicrosystems.us

German Facility:
LABO Microsystems GmbH
Bäckerstrasse 5
21244 Buchholz, Germany
Tel: + 49 4181 9906360
Fax: + 49 4181 9906369
E-mail: hq@labomicrosystems.com

Global Sales:
LABO Microsystems
P.O. Box 213671, Jumeirah, Dubai.
United Arab Emirates.
Tel: + 97 150 242 0048
Tel: + 91 941 666 6666
E-mail: labo@labomicrosystems.com



JAS-ANZ



ISO 9001:2000
File No. 18292.04

